

US Amendment to Technical Specification on page 26 and Symbols Used on page 27 in the Operating Manual (900751 071119)

IV Efficiency Level

Adaptor HemoCue HCA01 or Friwo FW7333/SM12 can also be used, input 100V~ - 240V~/50-60 Hz/200mA (<500mA, HCA01).



US Amendment to Technical Specification on page 26 and Symbols Used on page 27 in the Operating Manual (900751 071119)

IV Efficiency Level

Adaptor HemoCue HCA01 or Friwo FW7333/SM12 can also be used, input 100V~ - 240V~/50-60 Hz/200mA (<500mA, HCA01).

